matopat

Medisorb ABSORB

en

INSTRUCTIONS FOR USE
Medisorb ABSORB
Superabsorbent dressing non adherent

PRODUCT DESCRIPTION
Medisorb ABSORB superabsorbent dressing non adherent is an absorbent wound dressing that helps maintain a moist environment on the wound surface. It comprises an absorbent core layer that absorbs and retains exudate, which is wrapped between a hydrophobic breathable layer and a liquid permeable layer. The superabsorbent dressing absorbs and retains exudate and traps it within the absorbent core, thereby reducing the risk of maceration, and promoting wound healing. Medisorb ABSORB superabsorbent dressing non adherent has three layer structure which consists of: ➤ a hydrophobic breathable layer ➤ an absorbent core layer ➤ a liquid permeable layer.

Dressing material: polypropylene, polyethylene, cellulose, polyacrylate.

INTENDED PURPOSE
Medisorb ABSORB superabsorbent dressing non adherent is intended for use on moderate to heavily exuding wounds, such as: ➤ leg ulcers ➤ pressure ulcers ➤ diabetic foot ulcers ➤ surgical wounds.

CONTRAINDICATIONS
Do not use this product on dry or low exuding wounds, eyes, mucous membranes or in wound cavities because the dressing swells during absorption. Do not use it on third-degree burns. Medisorb ABSORB superabsorbent dressing non adherent should not be used on persons with known allergy or hypersensitivity to any of its components.

SIDE EFFECTS
Even though it may be rare, this product can have side effect, which may not occur for everybody. If any of the following side effects occur, stop using the product and consult your physician: ➤ Allergic (e.g. erythema, itching) reaction ➤ Pain ➤ Irritation

PRECAUTIONS AND WARNING
➤ The product is designed for single-using devices, in order to reuse them, may seriously damage their integrity and their performance, and cross-contamination may occur. ➤ Sterile. Do not use if sterile packaging is damaged or opened prior to use. Do not re-sterilize. ➤ Do not cut the dressing, the cover must remain intact. ➤ In the absence of available data supporting the use of this dressing in pregnant women, breast-feeding women, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used only after a risk-benefit assessment by the physician. ➤ Medisorb ABSORB superabsorbent dressing non adherent is not recommended to stand on the dressing, e.g. in the case of wounds on the feet. ➤ If the product does not begin to show signs of healing or if any other unexpected symptoms occur (e.g. fever or the wound or surrounding skin becoming red, warm or swollen), consult a health care professional for appropriate treatment.

PATIENT INTENDED POPULATION
The people who suffer leg ulcers, pressure ulcers, diabetic foot ulcers, surgical wounds.

MEDICAL CONDITION
Wound exudate

INTENDED USER OF THE DEVICE
Health care professionals, or by persons under supervision of health care professionals

INSTRUCTIONS FOR USE
Before using the dressing:
➤ Any excess hair should be clipped to ensure close approximation to the wound. ➤ Gently cleanse the wound and the surrounding skin in accordance with clinical practice. ➤ Dry the surrounding skin thoroughly.

Dressing selection and application
➤ Select a suitable size product to make sure the absorbent inner core overlaps the wound margins at least 1-2 cm. Do not trim the product. Do not fold the device. Folding the product may affect exudate uptake. ➤ Medisorb ABSORB superabsorbent dressing non adherent can be applied as the used both as a primary dressing or as a secondary dressing. When used as a primary dressing, Medisorb ABSORB superabsorbent dressing non adherent should be in close contact with the wound bed to ensure wound fluid is absorbed. ➤ 3. Secure the dressing after selecting an appropriate size. Medisorb ABSORB superabsorbent dressing non adherent can be secured in position with any other suitable fixation device, e.g. bandage, tape, net bandage. ➤ Dressing, the expansion during absorption into consideration. ➤ 4. Medisorb ABSORB superabsorbent dressing non adherent can be used in conjunction with compression therapy for venous leg ulcers.

Dressing change and removal
➤ 1. Medisorb ABSORB superabsorbent dressing non adherent may be left in place for several days depending on the condition of the wound and the surrounding skin, or as indicated by clinical practice. ➤ 2. Remove gently to avoid trauma to the wound bed. Remove in accordance with clinical practice.

STORAGE AND DISPOSAL
➤ 1. Store at room temperature in a place in specified temperature range and away from sunlight. ➤ 2. Dispose of in accordance with local environmental procedures.

SHELF LIFE – 3 years from date of manufacture.

OTHER INFORMATION
If any serious incident occurred in relation to the use of Medisorb ABSORB superabsorbent dressing non adherent, it should be reported to the manufacturer and competent authority of the Member State.

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