matopat

Medisorb ABSORB SILICONE

Dressing selection and application
• 1. Select a dressing size large enough to make sure the absorbent inner core overlaps the wound margins at least 1 cm. Do not cut the product. Do not fold the device. Folding the product may affect maximum absorbability. • 2. For dressings with three-dimensional shape, position the dressing so the middle release film first. Apply the adhesive side to the wound, allowing initial fixation of the dressing. For dressings with two release films, peel off the first release film and apply the adhesive side to the wound, allowing initial fixation of the dressing. • 3. Remove the second release film and smooth down the border on the skin. Do not stretch the dressing. • 4. Medisorb ABSORB SILICONE superabsorbent dressing adherent can be used both as a primary dressing or as a secondary dressing. • 5. Medisorb ABSORB SILICONE superabsorbent dressing adherent is not to be used in combination with compression therapy for various leg ulcers.

Dressing change and removal
• 1. Medisorb ABSORB SILICONE superabsorbent dressing adherent may be left in place for several days depending on the condition of the wound or surrounding skin, or as indicated by clinical practice. • 2. To remove the Medisorb ABSORB SILICONE superabsorbent dressing adherent, lift one corner and gently pull the dressing back on itself. • 3. Dressing removal can be facilitated by moistening the dressing with sterile saline.

STORAGE AND DISPOSAL
Store at room temperature (15-30°C) in specified temperature range and away from sunlight. • 2. Disposal should be handled according to local regulations.

SHELF LIFE – 3 years from date of manufacture.

OTHER INFORMATION
If any serious incident has occurred in relation to the use of Medisorb ABSORB SILICONE superabsorbent dressing adherent, it should be reported to the manufacturer and competent authority of the Member States.

en

INSTRUCTIONS FOR USE
Medisorb ABSORB SILICONE
Superabsorbent dressing adherent

PRODUCT DESCRIPTION
Medisorb ABSORB SILICONE superabsorbent dressing adherent is an absorbent wound dressing that helps maintain a moist environment on the wound surface. It comprises an absorbent core layer that absorbs and retains exudate, which is wrapped between a hydrophobic breathable layer and a liquid impermeable outer layer. The absorbent core contains superabsorbent fibers and particles. The Super Absorbent Dressing absorbs and retains exudate and **HARN** within the absorbent core, thereby reducing the risk of maceration, and promoting wound healing. The silicone wound contact layer facilitates the application of the dressing, and minimizes patient pain and trauma to the wound and surrounding skin when the dressing is removed.

INTENDED USE – Medisorb superabsorbent dressing adherent consists of: ▪ a hydrophobic breathable layer ▪ absorbent core ▪ liquid impermeable outer layer ▪ silicone wound contact layer. Dressing material contents: polypropylene, polyethylene, cellulose, polyacrylate, silicone.

INTENDED PURPOSE
Medisorb ABSORB SILICONE superabsorbent dressing adherent is intended for use on moderate to heavily exuding wounds, such as: ▪ leg ulcers ▪ pressure ulcers ▪ diabetic foot ulcers ▪ surgical wounds.

CONTRAINDICATIONS
Do not use the product on dry or low exuding wounds, eyes, mucous membranes or in wound cavities that require packing or wound filling absorption. Do not use if on third-degree burns. Medisorb ABSORB SILICONE superabsorbent dressing adherent should not be used on persons with known allergy or hypersensitivity to any of its components.

SIDE EFFECTS
Even though it may be rare, this product can have side effect, which may not occur for everybody. If any side effect occurs, stop using the product and consult your professional ▪ local reaction ▪ pain ▪ irritation ▪ itching ▪ allergic reaction

PRECAUTIONS AND WARNING
▪ The product is designed for single-processing devices, in order to reuse them, may seriously damage their integrity and their performance, and cross-contamination may occur. ▪ Sterile. Do not use if sterile packaging is damaged or opened prior to use. Do not re-sterilize. ▪ Do not use if the sterile barrier has been compromised. ▪ In the absence of available data supporting the use of this dressing on sensitive wounds, such as fistulas, the dressing should be used only after careful clinical assessment. ▪ Do not, on the contrary, on these applications. ▪ Do not use on wounds with visible signs of infection without medical recommendation. ▪ Do not use Medisorb ABSORB SILICONE – superabsorbent dressing adherent, in cavities. ▪ It is not recommended to stand on the dressing, e.g. in the case of wounds on the feet. ▪ If the wound does not begin to show signs of healing or if any other unexpected symptoms occur (e.g. signs of infection, excessive swelling) skin becoming red, warm or swollen), consult a health care professional for appropriate treatment.

INTENDED PATIENT POPULATION
The people who suffer leg ulcers, pressure ulcers, diabetic foot ulcers, surgical wounds.

MEDICAL CONDITIONS – Hospital and home setting.

INTENDED USER OF DEVICE
Health care professionals or lay persons under supervision of health care professionals.

INSTRUCTIONS FOR USE
Before using the dressing
• 1. Any excess hair should be clipped to ensure close approximation to the wound. • 2. Gently cleanse the wound and the surrounding skin in accordance with clinical practice. • 3. Dry the surrounding skin thoroughly.

Dressing selection and application
• 1. Select a dressing size large enough to make sure the absorbent inner core overlaps the wound margins at least 1 cm. Do not cut the product. Do not fold the device. Folding the product may affect maximum absorbability. • 2. For dressings with three-dimensional shape, position the dressing so the middle release film first. Apply the adhesive side to the wound, allowing initial fixation of the dressing. For dressings with two release films, peel off the first release film and apply the adhesive side to the wound, allowing initial fixation of the dressing. • 3. Remove the second release film and smooth down the border on the skin. Do not stretch the dressing. • 4. Medisorb ABSORB SILICONE superabsorbent dressing adherent can be used both as a primary dressing or as a secondary dressing. • 5. Medisorb ABSORB SILICONE superabsorbent dressing adherent is not to be used in combination with compression therapy for various leg ulcers.

Dressing change and removal
• 1. Medisorb ABSORB SILICONE superabsorbent dressing adherent may be left in place for several days depending on the condition of the wound or surrounding skin, or as indicated by clinical practice. • 2. To remove the Medisorb ABSORB SILICONE superabsorbent dressing adherent, lift one corner and gently pull the dressing back on itself. • 3. Dressing removal can be facilitated by moistening the dressing with sterile saline.

STORAGE AND DISPOSAL
Store at room temperature (15-30°C) in specified temperature range and away from sunlight. • 2. Disposal should be handled according to local regulations.

SHELF LIFE – 3 years from date of manufacture.

OTHER INFORMATION
If any serious incident has occurred in relation to the use of Medisorb ABSORB SILICONE superabsorbent dressing adherent, it should be reported to the manufacturer and competent authority of the Member States.

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