matopat

Medisorb P SILICONE BORDER

en

INSTRUCTIONS FOR USE
Medisorb P SILICONE BORDER
Soft silicone dressing

DESCRIPTION
Medisorb P SILICONE BORDER silicone foam dressing border is highly absorbent dressing that absorbs exudates and maintains a moist wound environment. The multi-layered construction facilitates dynamic fluid management to help to reduce the risk of maceration. The gentle silicone layer can be lifted and repositioned without losing its adherent. Also, this layer can minimize pain to patients and trauma to wounds and the surrounding skin. The waterproof outer layer protects the wound from dirt and bacteria.
The Medisorb P SILICONE BORDER silicone foam dressing border has an advanced five-layer structure which consists of: ► a soft silicone wound contact layer ► a polyurethane foam layer ► a non-woven distributing layer with technology ► a super absorbent polyacrylate fiber layer with technology ► a waterproof, breathable polyurethane outer layer.
Dressing material: silicone, polyacrylate, polyurethane, cellulose and polyolefin.

INTENDED PURPOSE
Medisorb P SILICONE BORDER silicone foam dressing border is designed for a wide range of exuding wounds such as: ► pressure ulcers ► leg and foot ulcers ► traumatic wounds ► surgical wounds ► skin tears.
The dressing may also be used as part of a prophylactic therapy to help prevent pressure ulcers.

CONTRAINDICATIONS
SIDE EFFECTS
Even though it may be rare, this product can have side effect, which may not occur for everybody. If any of the following side effects occur, consult a health care professional: ► local reaction (e.g. irritation, allergic reaction) ► pain ► infection

PRECAUTIONS FOR USE
► It should not be used on persons with a known allergy or hypersensitivity to the dressing or its components. ► Do not use if the sterile packaging is damaged or indicating agents such as hypochlorite solution or hydrogen peroxide. ► If the wound does not begin to show signs of healing or if any other unexpected reaction occurs, stop the use of the dressing and consult a health care professional. ► Consult with a health care professional about appropriate compression bandages. ► Do not use on wounds with exposed bone, tendon or muscle. ► Do not use on wounds with extensive contamination may occur. ► Sterile. Do not use if sterile packaging is damaged or opened. ► Single use only. ► Do not use after the expiration date. ► Always use dressings supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, elderly and people at risk. ► Refer to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.

INTENDED PATIENT POPULATION
The people who suffer pressure ulcers, leg and foot ulcers, traumatic wounds, surgical wounds, skin tears.

MEDICAL CONDITIONS – Hospital and home setting.

INTENDED USER OF THE DEVICE
Health care professionals or any person under supervision of health care professionals.

INSTRUCTIONS FOR USE
Before using the dressing:
► 1. Any excess hair should be clipped to ensure close approximation to the wound. ► 2. Gently cleanse the wound and surrounding skin with normal saline solution or with clinical practice. ► 3. Dry the surrounding skin thoroughly.
Dressing selection and application:
► 1. Select an appropriate absorptive. For difficult to dress anatomical locations, such as the sacrum,

heel or multisite, the specifically shaped adhesive dressing may be used. For use on wounds: the wound pad should be dry to accommodate exudate. ► 2. Peel off the backing paper and apply the dressing lightly and smooth over the border so that the skin Oz not stretch the dressing. ► 4. When necessary, fixate the dressing with bandage, tapes, tubular bandage, straps and the use of appropriate and commonly used clinical practice.
Dressing change and removal:
► 1. The recommended change of Medisorb P SILICONE BORDER silicone foam dressing border may be left in place for several days depending on the condition of the wound and the surrounding skin, or as indicated by clinical practice. For use on prevention: the area at risk of pressure ulcerSusue damage should be inspected at regular intervals and dressing changed as required. ► 2. Remove Medisorb P SILICONE BORDER silicone foam dressing border, lift one corner and gently pull the dressing back on itself.

STORAGE AND SHELF LIFE
► 1. Keep in a clean and dry place in specified temperature range and away from sunlight. ► 2. Dispose product should be handled according to environmental procedures.

SHELF LIFE – 3 years from date of manufacture.

OTHER INFORMATION
If any serious incident has occurred in relation to the use of Medisorb P SILICONE BORDER silicone foam dressing border, it should be reported to the manufacturer and competent authority of the Member States.

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